Stability of Antisera and Reagents during Shipment

Our antisera and reagents are shipped at room temperature without cooling.

The quality of the products is not affected during transport, even at elevated summer temperatures.

After receipt, the products should be kept at the indicated temperature for long-term storage:

IF-antisera and secondary antibodies: at or below -20°C (avoid repeated freezing- and thawing cycles, freeze in suitable aliquots, if applicable)

ELISA reagents and controls: +4°C (refrigerator)

PCR reagents: at or below -20°C

Buffers and chemicals: as indicated on the product labels

LOEWE® Biochemica GmbH